



Product Datasheet

ATP6V1C2 antibody (orb400680)

Description	ATP6V1C2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human V-type proton ATPase subunit C 2 protein (190-262AA)
Target	ATP6V1C2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NEY4
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat heart tissue...



Immunofluorescence analysis of A549 cell...