



Product Datasheet

ACY3 antibody (orb400700)

Description	ACY3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human N-acyl-aromatic-L-amino acid amidohydrolase (carboxylate-forming) protein (5-319AA)
Target	ACY3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96HD9
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cel...



Immunofluorescence analysis of HepG2 cel...