

Product Datasheet

Pdk3 antibody (FITC) (orb400722)

Description	Pdk3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 3, mitochondrial protein (131-362AA)
Target	Pdk3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q922H2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.