



Product Datasheet

HSD3B7 antibody (orb400732)

Description	HSD3B7 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human 3 beta-hydroxysteroid dehydrogenase type 7 protein (131-197AA)
Target	HSD3B7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H2F3
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver tis...



Immunofluorescence analysis of HepG2 cel...