



Product Datasheet

ACP6 antibody (orb400756)


Description	ACP6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Lysophosphatidic acid phosphatase type 6 protein (35-324AA)
Target	ACP6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NPH0
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of PC3 whole cell ...



Immunohistochemical staining of human pr...



Immunohistochemical staining of Human ki...