




Product Datasheet

TBX1 antibody (orb40078)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TBX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43133
Uniprot ID	O43435
NCBI	NP_005983.1 , NP_542378.1 , NP_542377.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of paraffin-embedded tissue.