

Product Datasheet

CD55 antibody (FITC) (orb400818)

Description	CD55 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Complement decay-accelerating factor protein (35-126AA)
Target	CD55
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08174
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.