

Product Datasheet

ND2 antibody (orb40082)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target MT-ND2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG



Clonality Polyclonal

MW 38961


Uniprot ID [P03891](#)

Dilution Range WB 1:50-1:100

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of ND2 antibody...