



Product Datasheet

TNFRSF10A antibody (orb400932)


Description	TNFRSF10A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 10A protein (263-468AA)
Target	TNFRSF10A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000220
Dilution Range	IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunohistochemical staining of Human ki...



Immunohistochemical staining of human CO...