

Product Datasheet

bmp2b antibody (Biotin) (orb400963)

Description	bmp2b antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Bone morphogenetic protein 2 protein (24-411AA)
Target	bmp2b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	093369
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.