

Product Datasheet

tal1 antibody (HRP) (orb400965)

Description	tal1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish T-cell acute lymphocytic leukemia protein 1 homolog protein (1-324AA)
Target	tal1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O93507
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.