


Product Datasheet

btFP antibody (orb400988)

Description	btFP antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bacteroides fragilis Fragilysin protein (26-405AA)
Target	btFP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54355
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...