


Product Datasheet

phr antibody (orb400992)

Description	phr antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Synechococcus sp. Deoxyribodipyrimidine photo-lyase protein (2-484AA)
Target	phr
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05327
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant phr...