


Product Datasheet

THI2.1 antibody (orb400996)

Description	THI2.1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Viscum album Viscotoxin-A3 protein (27-72AA)
Target	THI2.1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01538
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...