

Product Datasheet

ghrA antibody (FITC) (orb401006)

Description	ghrA antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Glyoxylate/hydroxypyruvate reductase A protein (1-312AA)
Target	ghrA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7LFE3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.