



Product Datasheet

PDK1 antibody (orb401032)

Description	PDK1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial protein (294-432AA)
Target	PDK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15118
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse kidney ti...



Immunohistochemical staining of human li...