


Product Datasheet

TAF8 antibody (orb401056)

Description	TAF8 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Transcription initiation factor TFIID subunit 8 protein (122-311AA)
Target	TAF8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z7C8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat kidney tiss...