


Product Datasheet

Fuca1 antibody (orb401076)

Description	Fuca1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Tissue alpha-L-fucosidase protein (18-452AA)
Target	Fuca1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99LJ1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat kidney tissue...