

Product Datasheet

SENP2 antibody (HRP) (orb401121)

Description	SENP2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Sentrin-specific protease 2 protein (138-324AA)
Target	SENP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HC62
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.