

Product Datasheet

Snaclec rhodocytin subunit alpha antibody (HRP) (orb401125)

Description	Snaclec rhodocytin subunit alpha antibody (HRP)
Species/Host	Rabbit
Reactivity	Reptile
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Calloselasma rhodostoma Snaclec rhodocytin subunit alpha protein (1-136AA)
Target	Snaclec rhodocytin subunit alpha
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9I841
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.