



Product Datasheet

PSD antibody (orb401228)


Description	PSD antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human PH and SEC7 domain-containing protein 1 protein (247-489AA)
Target	PSD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A5PKW4
Dilution Range	WB:1:500-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cel...



Immunohistochemical staining of human br...



Immunohistochemical staining of human sm...