

Product Datasheet

HLA-C antibody (FITC) (orb401282)

Description	HLA-C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human HLA class I histocompatibility antigen, Cw-7 alpha chain protein (62-306AA)
Target	HLA-C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10321
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.