



Product Datasheet

ASPM antibody (orb401324)


Description	ASPM antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Abnormal spindle-like microcephaly-associated protein (101-234AA)
Target	ASPM
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IZT6
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal gland tissue showing ASPM expression.



Immunohistochemical staining of human lung tissue showing ASPM expression.



Immunofluorescence analysis of Hela cells showing ASPM expression.