



Product Datasheet

NMT2 antibody (orb40133)


Description	NMT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Recombinant Human Glycylpeptide N-tetradecanoyltransferase 2 protein (1-120AA)
Target	NMT2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60551
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of HepG2 wh...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human te...