



Product Datasheet

MARS2 antibody (orb401336)


Description	MARS2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Methionine--tRNA ligase, mitochondrial protein (225-460AA)
Target	MARS2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96GW9
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human br...