

Product Datasheet

MCU antibody (FITC) (orb401358)

Description	MCU antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Calcium uniporter protein, mitochondrial protein (155-233AA)
Target	MCU
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NE86
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.