

Product Datasheet

TAB3 antibody (HRP) (orb401361)

Description	TAB3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human TGF-beta-activated kinase 1 and MAP3K7-binding protein 3 protein (361-712AA)
Target	TAB3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N5C8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.