



Product Datasheet

SEMA5B antibody (orb401376)


Description	SEMA5B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Semaphorin-5B protein (845-1013AA)
Target	SEMA5B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P283
Dilution Range	IHC-P:1:20-200, IF:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart tis...



Immunohistochemical staining of human te...



Immunohistochemical staining of human lu...