



Product Datasheet

PACSIN3 antibody (orb401408)


Description	PACSIN3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Protein kinase C and casein kinase substrate in neurons protein 3 protein (1-300AA)
Target	PACSIN3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKS6
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis
of Rat liver tissu...



Immunohistochemical
staining of human
ad...



Immunohistochemical
staining of human
br...