

Product Datasheet

notch1a antibody (orb401428)

Description	notch1a antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Neurogenic locus notch homolog protein 1 protein (21-338AA)
Target	notch1a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46530
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.