

Product Datasheet

tbx5a antibody (orb401436)

Description	tbx5a antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish T-box transcription factor TBX5-A protein (1-492AA)
Target	tbx5a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9IAK8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.