



Product Datasheet

MAP3K9 antibody (orb401516)

Description	MAP3K9 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 9 protein (864-1098AA)
Target	MAP3K9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P80192
Dilution Range	IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat lung tissue...



Immunohistochemical staining of human lung tissue...