


Product Datasheet

MNX1 antibody (orb401524)

Description	MNX1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Motor neuron and pancreas homeobox protein 1 protein (50-99AA)
Target	MNX1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50219
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...