


Product Datasheet

OLIG1 antibody (orb401536)

Description	OLIG1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Oligodendrocyte transcription factor 1 protein (17-105AA)
Target	OLIG1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TAK6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat heart tissue...