

Product Datasheet

OPN1SW antibody (FITC) (orb401546)

Description	OPN1SW antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Short-wave-sensitive opsin 1 protein (307-348AA)
Target	OPN1SW
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03999
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.