


Product Datasheet

PTGER4 antibody (orb401560)

Description	PTGER4 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Prostaglandin E2 receptor EP4 subtype protein (354-472AA)
Target	PTGER4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35408
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat liver tissue...