


Product Datasheet

Syntaxin 1A antibody (orb401572)

Description	Syntaxin 1A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Syntaxin-1A protein (1-265AA)
Target	STX1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16623
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of PC3 whole cell ...