



Product Datasheet

PDE4B antibody (orb401624)

Description	PDE4B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human cAMP-specific 3',5'-cyclic phosphodiesterase 4B protein (27-215AA)
Target	PDE4B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07343
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human glioma tissue.