

Product Datasheet

ROR2 antibody (Biotin) (orb401631)

Description	ROR2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tyrosine-protein kinase transmembrane receptor ROR2 protein (752-926AA)
Target	ROR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01974
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.