

Product Datasheet

HLA-C antibody (FITC) (orb401646)

Description	HLA-C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human HLA class I histocompatibility antigen, Cw-3 alpha chain protein (200-306AA)
Target	HLA-C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04222
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.