

Product Datasheet

pbpA antibody (Biotin) (orb401651)

Description	pbpA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Clostridium botulinum Penicillin-binding protein 1A protein (663-830AA)
Target	pbpA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A5I6G4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.