

Product Datasheet

PAG 1 antibody (FITC) (orb401682)

Description	PAG 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	FITC
Immunogen	Recombinant Bovine Pregnancy-associated glycoprotein 1 protein (54-380AA)
Target	PAG 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q29432
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.