



Product Datasheet

FZD4 antibody (orb40172)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FZD4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59881
Uniprot ID	Q9ULV1
NCBI	NP_036325.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in mouse heart

...