



Product Datasheet

RAET1E antibody (orb401732)


Description	RAET1E antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Retinoic acid early transcript 1E protein (31-225AA)
Target	RAET1E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3ZAR2
Dilution Range	IHC-P:1:20-200, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...



Immunohistochemical staining of human sk...



Immunofluorescence analysis of HepG2 cel...