

Product Datasheet

WWOX antibody (HRP) (orb401749)

Description	WWOX antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human WW domain-containing oxidoreductase protein (1-180AA)
Target	WWOX
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZC7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.