

Product Datasheet

RPS10P5 antibody (FITC) (orb401754)

Description	RPS10P5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Putative 40S ribosomal protein S10-like protein (1-176A)
Target	RPS10P5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NQ39
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.