



Product Datasheet

GFAP antibody (orb401780)


Description	GFAP antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human GFAP protein (292-432AA)
Target	GFAP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P14136
Dilution Range	WB: 1:500:5000, IHC-P: 1:50:500, IF/ICC: 1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SY5Y cell...



Immunofluorescence analysis of U251 cell...



Immunofluorescence analysis of U87 cells...