



## Product Datasheet

### BFAR antibody (orb40184)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | BFAR antibody                                                        |
| <b>Species/Host</b>        | Rabbit                                                               |
| <b>Reactivity</b>          | Human                                                                |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Tested Applications</b> | ELISA, IHC                                                           |
| <b>Immunogen</b>           | Recombinant Human Bifunctional apoptosis regulator protein (1-140AA) |
| <b>Target</b>              | BFAR                                                                 |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                    |
| <b>Form/Appearance</b>     | Liquid                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.        |
| <b>Note</b>                | For research use only                                                |
| <b>Isotype</b>             | IgG                                                                  |
| <b>Clonality</b>           | Polyclonal                                                           |
| <b>Uniprot ID</b>          | <a href="#">Q9NZS9</a>                                               |
| <b>Dilution Range</b>      | IHC-P:1:20-200                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                      |



Immunohistochemical staining of human sm...



Immunohistochemical staining of human br...