

Product Datasheet

SF2523 (orb401842)

Description	SF2523 is a highly selective and potent inhibitor of PI3K and BRD4 with IC50 values of 16 and 241 nM ...
Note	For research use only
Formula	$C_{19}H_{17}NO_6S$
MW	371.41
CAS Number	[1174428-47-7]
SMILES	<chem>C1COCCN1C2=CC(=O)C3=C(O2)C(=CS3)C4=CC5=C(C=C4)OCCO5</chem>
Expiration Date	12 months from date of receipt.