

Product Datasheet

TMP195 (orb401905)

Description	TMP195 is the most potent and selective class IIa HDAC inhibitor.
Note	For research use only
Formula	$C_{32}H_{19}F_3N_4O_3$
MW	456.42
CAS Number	[1314891-22-9]
SMILES	<chem>CC(C)(CNC(=O)C1=CC=CC(=C1)C2=NOC(=N2)C(F)(F)F)C3=COC(=N3)C4=CC=CC=C4</chem>
Expiration Date	12 months from date of receipt.